Mathis

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Application

Textiles, nonwovens, foils, papers etc. are coated mostly by spreading a coating paste regularly onto the carrying material by means of a doctor knife. In order to enable laboratory tests our LABCOATER has been developed, which mainly consists of the following integral parts:

1. Dryer
2. Guiding and supporting frame
3. Coating head
4. Pin frame (material holder)

By means of this coating combination it is possible to introduce and dry the coated goods directly in the dryer without re-pinning. Therefore, also multiple coatings can be realised quickly and reliably.

Range of Application:

The coating device enables following coating methods:

- Floating knife
- Roller knife
- Rubber blanket knife.

Beside knife application other accessories enable also:

- Lamination

1. Dryer

This unit mainly consists of following parts:

Sample preparation

The samples of max. 33×43 cm are fixed on a pin frame, adjustable in length or in width and length, whereby the fabric tension can be varied between 0–14 kp. The pin frame transport is functioning automatically. In both end positions, in and out, the entrance slit is closed.

Process data

The required entries for the process details such as dwell time, temperature and fan speed are entered at the touch screen of the process controller.

Concept

An efficient circulation fan, infinitely variable, allows selecting optimum treatment conditions for various materials. The use of a suitable air channel and corresponding outlet surfaces enable an even ventilation on both sides.

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Temperature measurement and heating

The temperature measurement is based on the resistance measurement by means of Pt 100. In order to have a precise value, the temperature sensor is placed shortly before the air outlet.

The heating elements are located above the circulation fan at the re-entry into the air channel. In this way the circulation air is reheated shortly before the re-entry into the air channel.

By an additional temperature sensor, the maximum temperature of the heating elements is measured and also limited to a maximum value of 260 °C in order to protect the heaters.

Small space requirement

By the use of the process control unit, several operating elements can be omitted and therefore the corresponding space requirement is considerably reduced. Nevertheless the electric components are well accessibly placed underneath of the drying chamber. Easily removable side panels allow a fast access if necessary.

Transport frame

The entry of the material holder is activated by pressing a corresponding button and after running down of the selected dwell time the sample is delivered automatically.

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Reliability

Several supervising systems are installed into this unit in order to increase the reliability, for instance the above mentioned temperature control for the heating. Furthermore the fan is equipped with a rotation control which cuts off the heating in case of problems with the fan motor or breakage of the fan belt.

Accessibility to the inside of the dryer

A good accessibility is given by the removable cover, which can be lifted upwards.

Heat insulation

A good heat insulation is obtained by a double wall construction with a well dimensioned insulation thickness.

Cooling

A quick cooling can be reached by an easy lifting of the cover. On request the apparatus can be delivered with a lifting motor for the cover.

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Options for the dryer

- **IR radiation pyrometer**
For the contactless measurement of the sample surface temperature. Several different treatments can be made: Pre setting of an air temperature and measuring of the surface temperature or required treatment temperature starts only when the set surface temperature is reached.
- **Exhaust fan**
An exhaust fan is available for heat treatments which cause some fumes or even heavy smoke, and therefore require such a unit. (option)
- **Positioning device for air guiding flap**
On the standard version the position of the air flap for the distribution of the top/bottom air flow is made by a star shaped handle and the position is fixed with a knurled screw. As an option the flap position can be varied by an electric motor. Advantage: better air and temperature distribution as well as easier reproduction of settings.

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Process control

The programmable process controller Mathis UNIVISION Touch S is very user friendly and simple to use, since a lot of information is presented in text form. The entry of process data such as dwell time, temperature, fan speed etc. is done via menu selection with English or German text at the touch screen panel. Beside entries for a single treatment consisting of temperature, fan speed and following dwell time, also an automatic program sequence consisting of several subsequent treatments at different temperatures and dwell times can be programmed.

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Daily or weekly program

Apart from entries for several subsequent treatments which are automatically executed, also complete daily or even full weekly programs can be entered at the programmer. In such kind of working plan also stand-by modus can be entered. In case of using the apparatus in double or multiple shift operation, the program can be set that the unit is ready for work at the required process conditions such as temperature, etc. set for Monday early morning.

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Display

On a clearly visible display the entered process data such as dwell time, temperature, fan speed as well as status information of the functions such as heating, cooling are permanently visible. After entry of the sample these values are permanently updated. In case of operating errors or possible machine errors, a corresponding message will be presented. A comprehensive help file is also available.

Set up of process data is done either in dialog mode with the touch screen or can be done at a personal computer with corresponding software.

Entries

All the process data, machine parameters as well as other information are saved onto a MMC memory card. By the use of this card the data can be transferred to other controllers of the same type. The same controller enables also a great number of inter lockings of various machine functions. This reduces operating errors to a great extent. Also a great number of function controls can be realised for example for temperature sensor control etc. which enables a good reliability. In case of errors which might happen during the process an alarm signal as well as a corresponding error statement will be presented at the display. This text information also includes some statements about possible remedy concerning the actual error.

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
Process control via personal computer

If this option is required the data can be transferred to a personal computer by the MMC memory card. There is also a possibility of a data exchange of the controller directly to the PC via Interface RS 232 or by means of a Ethernet interface for TCP/IP connection with a LAN network. This enables to check the process conditions of the actual procedure.

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
Protocol printer

A protocol print out for example for ISO 9001 is possible in one of the above mentioned connections with a PC and its linked printer. The print out can be used for quality control of the actually processed data including possible error statements. The attachment of the frame including coating device can also be made later on.

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Example for temperature recording of air temperature and IR surface temperature

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2. Guiding and supporting frame

This frame serves also as a guiding device for the coating head and is attached to the lower part of the dryer.

The attachment of this frame including coating device can also be made later on.

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3. Coating Head

For coating, the guiding and supporting frame is lifted that the supporting rollers touch the carrying material. In this way the coating head is in working position.

Depending on the rollers or rubber blanket employed, it is possible to work either with floating knife, roller or rubber blanket knife.

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The doctor knife is added vertically, the angle position being adjustable. The coating thickness can be adjusted and read off from clock gauges (accuracy 0.01 mm). By means of two handles the coating head is guided over the carrying material. At the end of the spreading procedure, the coating head is to lower again so that the pin frame can enter into the dryer.

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Options for coating head

- **Motor drive for coating head**
The forward movement of the coating head which is made by pulling at the two handles can on request also be made with a motor drive. The coating speed of the forward movement is infinitely variable.
- **Different knife shapes**
Suitable coatings knives are available for PVC and PU coatings.

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


Beside the above mentioned roller knife also floating knife application can be made with the standard accessories for the coating head.

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Additionally also rubber blanket knife application can be done by the use of an endless rubber blanket which is placed around the two floating rollers. (without picture)

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This additional part consists of a holding device with pressing roller plus material carrier. This device is placed at the coating head instead of the doctor knife.

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Application: lamination of carriers after coating of the adhesives e.g. PVC, dispersions etc. in the reversing method.

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Essentially the same as type SV-K however with fast lifting device of the laminating roller. Laminating is done against a separate roller. The complete fast laminating device is to be placed behind the doctor knife and fastened by means of a clamping lever. This enables a laminating of the carrier immediately after spreading of the adhesive.

Application: laminating of carriers immediately after coating of the adhesive, especially solvent containing adhesives.

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To stretch material, which can be pinned up. The tension is adjustable in length only. And can be read off on a scale.

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Pin frame adjustable in length and in width

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To stretch material which can be pinned up. The tension is adjustable in length and in width. The pin frames are also available with longer needles, which make it easier to hold carpet samples.

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
Paper and foil frame

To hold release paper, foils non-woven fabrics etc. This can be done by means of a clamping device. The sample is tensioned in length.

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Reversing pin frame

On this special pin frame the material is tensioned in length. The inner part of the frame with the pinned material can be turned. In this way it is possible to coat both sides of the sample without repinning.

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Sieve frame

The sieve frame made of stainless steel wire netting with removable cover is specially suited to place loose material.

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Our manufacturing programme includes in addition:

Laboratory Padders vertical, horizontal
Thermolising range
Drying and curing range
Pad-steam range
Beaker-Dyer
Jet-Dyer
Bobbin-Dyer
Coating device
Jig
Apparatus for discontinuous coating
Continuous Coating ranges

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Werner Mathis AG
Werner Mathis U.S.A. Inc.
Concord N.C. 28026, U.S.A.