*Mathis
Thermotester
Typ LTE-T*

*Mathis
Mathis
Mathis
Mathis
Mathis*

*Thermotester
Type LTE-T*

Allgemein

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General

This thermotester is a largely automatic oven for the determination of the thermal (static) stability of plastics.

In combination with the laboratory coater and other accessories the thermotester can be used as a gelling oven as well as an universal testing oven in the application technique of plastics.

1. Thermotester LTE-T

The development of this apparatus was based upon the following points:

- 1. Fast and labour-saving preparation of the sample stripes and the apparatus.*
- 2. Accurate following of the chosen testing conditions in order to receive reproducible results for production machines.*
- 3. The most possible number of automatic operations concerning the test possibilities in order to keep the labour as small as possible.*

Furthermore it is especially remarkable that it is not necessary to open the oven during a test process. The sample stripes are automatically delivered through a narrow slit and are intensively cooled outside the thermic field. By the special design of the cooling zone the entrance of cool air into the oven is impossible.

Another important detail is also that sample stripes with thermal overload can be taken away from the apparatus without interruption of the running process in order to make sure that the other samples are not influenced by some decomposing products.

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Carrying frame incl. 14 U-shaped holders. All sample stripes have the same size and are placed onto the holders.

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Sieve part incl. holding clamps for sample stripes having all the same size or different sizes up to the full sieve area.

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1.1. Procedure of a thermo test

Heat the thermotester up to required temperature.

Place the 14 sample stripes onto the U-shaped sample holders. The latter are arranged in a holding frame side by side.

After entering into the thermotester the sample stripes remain in their end position at the adjusted temperature until the ground time has run down.

Then, the sample stripes (all 14 of them) are moved out by one interval distance. Now the interval time starts running. Each time the interval time has run down the sample stripes are moved out the same interval distance until the holding frame is completely out.

Now a line or zone can be recognized from where a darker colour indicates that the material is thermically overloaded or burned. Based on the selected ground time plus the x-times repeated interval distance (which represents each time the selected interval time) the total dwell time can be calculated which the tested material can remain under the adjusted temperature until a change in colour indicates a thermal overloading or burning.

1.2 Methods of testing

The thermotests can be made in following ways:

1. In order to get directly comparable results of different materials.

For this purpose several different materials are tested simultaneously.

2. In order to get finally determined values of the same material.

For this purpose a number of sample stripes of the same material are tested.

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
Great advantages are also given by the use of standardized ground time, interval distance and interval time.

This allows a direct comparison of new materials or stabilizers with formerly known test results of common, existing materials.

1.3 Basic concept

Central control unit

Instead of a control box this apparatus has a compact unit on which all the control parameters can be entered.

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These include:

- temperature
- ground time
- interval time
- interval distance
- fan revolutions

The necessary electric elements like transport motor, contactors, relays etc. are built in underneath of the drying chamber.

In order to guarantee a good heat insulation the main box is an insulated double-wall construction.

Temperature is controlled by an electronic measuring and regulating instrument within the range of 20–250 °C.

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The electric heating capacity is approximately 7.2 kW.

The circulation fan is driven by a 0.75 kW AC motor.

A quick cooling can be reached by an easy lifting of the cover. On request a lifting motor for the covering hood will be supplied.

The interior of the thermotester is accessible by simply taking off the cover.

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2. Thermotester with coating device (accessory) LTE-TS

On request the thermotester can be combined with a coating device which enables following applications:

- floating knife
- roller knife
- rubber blanket knife

With suitable accessories also following processes are possible:

- lamination
- roller application

2.1 Coating Head

For coating, the guiding and supporting frame is lifted that the supporting rollers touch the carrying material. In this way the coating head is in working position.

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The doctor knife is added vertically, the angle position being adjustable. The coating thickness can be adjusted and read off from clock gauges (accuracy 0.01 mm). By means of two handles the coating head is guided over the carrying material. At the end of the spreading procedure, the coating head is to lower again so that the pin frame can enter into the dryer.

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Motor drive for coating head

The forward movement of the coating head which is made by pulling on two handles can on request also be made with a motor drive. The moving speed is infinitely variable.

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

Different knife shapes

Suitable coating knives are available for PVC and PU coatings.

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2.2 Material holders

Pin frame adjustable in length

To stretch material which can be pinned up. The tension is adjustable only in length and can be read off on a scale.

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Pin frame adjustable in length and in width

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To stretch material which can be pinned up. The tension is adjustable in length and in width.

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The pin frames are also available with longer needles, which make it easier to hold carpet samples.

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
To hold release paper, foils, non-woven fabrics etc. This can be done by means of a clamping device. The sample is tensioned in length.

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
On this special pin frame the material is tensioned in length. The inner part of the frame with the pinned material can be turned. In this way it is possible to coat both sides of the sample without re-pinning.

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The sieve frame made of stainless steel wire netting with removable cover is specially suited to place loose material.

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2.3 Accessories

Laminating device type SV-K

This additional part consists of a holding device with pressing roller plus material carrier. This device is placed at the coating head instead of the doctor blade.

Application: laminating of carriers after coating of the adhesives e.g. PVC, dispersions etc. in the reversing method.

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Fast laminating device type K-SK

Essentially the same as type SV-K however with a fast lifting device of the laminating roller. Laminating is done against a separate roller. The complete fast laminating device is to be placed behind the doctor knife and fastened by means of a clamping lever. This enables a laminating of the carrier immediately after spreading of the adhesive.

Application: laminating of carriers immediately after coating of the adhesive, especially solvent containing adhesives.

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
Roller coating device type SV-WB


For carpet back coatings, coatings of velvet etc. by means of roller application. The coating roller is mechanically driven. There are two doctor blades, one for the application roller and the second one for the carpet.

The roller coating device is simply attachable. It is placed onto the holding device instead of the supporting rollers.

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* Zusatz

Thermotester LTE-T

| | |
|--------------------------------|-----------------------------|
| temperature range | 20–250 °C |
| heating | 7.2 kW |
| heating time | to 180 °C approx. 25 min |
| fresh and exhaust air | adjustable |
| entry of samples | automatic |
| no. of sample stripes | 14 |
| size of sample stripes | 2 x 50 cm |
| delivery distance per interval | 1–100 mm |
| interval time | 3 s–999 min |
| ground time | 3 s–999 min |
| weight | approx. 220 kg |

Thermotester with coating device LTE-TS

| | |
|------------------------|-----------------|
| size of sample | max. 33 x 43 cm |
| size of coated surface | max. 28 x 30 cm |
| coating speed* | 0.2–5 m/min |
| weight | approx. 240 kg |
| width | approx. 95 cm |

*Option

Mathis

Textilmaschinen Textile Machines
Laborapparate Laboratory Equipment
Beschichtungsanlagen Coating Ranges

Werner Mathis AG
Rütisbergstrasse 3
CH-8156 Oberhasli / Zürich
Switzerland

Telefon 41(0)1-852 5050
Telefax 41(0)1-850 6707
E-mail info@mathisag.com
Homepage www.mathisag.com

Tochterfirmen / sister companies:

Werner Mathis U.S.A. Inc.
2260 HWY 49 N.E. / P.O. Box 1626
Concord N.C. 28026, U.S.A.
Phone: 1-704-786-6157
Fax: 1-704-786-6159

Aparelhos de laboratorio Mathis Ltda
Estrada dos Estudantes 261
CEP 06700 COTIA S.P., Brasil
Phone: 55 (0)11-7922-2333
Fax: 55 (0)11-7922-2598
E-mail mathbra@ibm.net