CT3[™] Texture Analyzer

compression and tension testing for rapid QC analysis

An extensive history and customer input have contributed to the development of the most powerful, low cost, stand-alone Texture Analyzer ever produced. With six test modes (plus calibration check) and a wide choice of accessories, no other texture analyzer has ever done so much without requiring a computer and software!

Standard Test Modes

Normal Test: a single compression cycle Hold Time Test: compress and hold **Cycle Count Test:** compress multiple times **Bloom Test:** gelatin bloom strength test **TPA Test:** texture profile analysis **Tension Test:** tensile testing Surimi Test: gel strength Static Load Test: calibration check **Texture Loader Software**

allows up to ten custom tests and ability to lock parameters

Compression distance

up to 10cm, can accommodate sample up to 22.5cm, almost 9 inches tall. Probe shaft is 8cm from back wall.

Choice of Load Cells 7 measurement ranges up to 50kg

Choice of Base Tables allows for larger samples and more accessory choices



CT3 with Fixture Base Table and Cylindrical Probe in compression (TPA) mode

What's Included?

Instrument with choice of load cell Texture Loader Software USB & Power Cables

What else do I need?

Rotary Base Table, Fixture Base Table (see below) or Adjustable Base Table (p66) At least one probe or test fixture (p60-64)

Optional Accessories

The CT3 has a wide variety of probes, fixtures and jigs which enable it to be very versatile. AMETEK Brookfield can also custom design a fixture and probe for most applications.

TexturePro CT Software TA-CT-PRO-AY (p59) Temperature Probe DVP-94Y Bubble Level TA-LVL Calibration Weight Set (p59) Gelatin Bath System for gel conditioning (p57) Bloom Jar - industry approved TA-GBB-2

TexturePro CT Software Optional


COLLECT DATA AND PERFORM DETAILED DATA ANALYSIS WITH REAL-TIME GRAPHIC PLOTTING.

While the CT3 can perform many tests in stand alone mode, use of TexturePro CT Software permits creation of multiple tests and automatic execution without operator involvement.

Easily create custom reports and graphs right from the menu screen.

- Sample identification set-up screen helps new operators quickly get started; test fields outline a variety of parameters
- Intuitive set up for test methods and database file structures in a single window
- Data is captured as a graph and stored in tabular database format
- Advanced data analysis with built-in parameter calculations such as springiness, chewiness, hardness and much more!

thi Copp Loads RefutS C X B 20 Venteent Martial and A 20 Venteent Martial and A 20 Venteent Page 10 (place Page 10 (place Page 10 (place) Calcelar makine C Scolar and Martial Colar C Scolar Colardiane C Martial Colardiane C Val. C Val. C Val.	First Cycle Calculations □ Nardress □ N	Lotimutos a Tape Apparent Hobbers Apparent Hobbers Cartopose Cartopose Cartopose PableStess Pabless Adversement Adversement Paelence Tompers Unruph	Conset Cycle Calculation device and Cycle Calculation device and Cycle Calculation Seconstation Consets Seconstation Consets
	Querity of Fischeres Fracturability To: Fracture Load Drop 01	IntFacture Volk Done □ IntFacture Defumation □ IntFacture N Defumation	



Sample Test Set-up

MODEL	Load Range / Resolution*
CT3-100	0-100g/0.01g
CT3-1000	0-1000g/0.10g
CT3-1500	0-1500g/0.20g
CT3-4500	0-4500g/0.50g
CT3-10kg	1-10000g/1.0g
CT3-25kg	1-25000g/2.0g
CT3-50kg	2-50000g/5.0g

g = grams kg = kilograms *Accuracy = $\pm 0.5\%$ Full Scale Range (FSR)

On Screen Live Force Deformation Curve

ALL CT3 MODEL SPECIFICATIONS		
Speed:		
Range	0.01-0.1mm/s (increments 0.01mm/s)	
	0.1 - 10mm/s (increments 0.1mm/s)	
Accuracy	±0.1% of setspeed	
Position:		
Range	0-101.6mm	
Resolution	0.1mm*	
Accuracy	0.1mm	
mm = millimeter s = seconds		

*Resolution 0.01mm when used with TexturePro CT Software



Tension Mode

provides tensile testing capability



TA-CW-1500C

Calibration Weight Set contains a combination of certified weights which may be used to confirm the calibration and linearity of each specific load cell.