

基本型離心機 Basic Type



C-12B



掀蓋內部 Inside

特點：

- 適用於實驗室及臨床實驗。
- 可調速、操作容易。
- 6ml/10ml之試管可供選擇。

Feature：

- Suitable for laboratory and clinical tests
- Speed adjustable design, easy operations
- 6ml and 10ml adapter are optional

Model	C-12B
轉速 Speed (rpm)	0~2,500
最大離心力 Max. RCF (xg)	580 (6ml x 4pcs)
	650 (10ml x 4pcs)
最大容量 Max. Capacity (ml)	40
定時器 Timer (min:sec)	----
重量 Weight (kg)	3
尺寸 Size (W x L x H mm)	250 x 250 x 250
電源 Power	AC 110V / 220V, 50/60Hz

◎ 離心試管另購。 Tubes not included.

■ Rotor sets for C-12B



RA-0604



RA-1004

Rotor	Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
RA-0604	6ml x 4	10.5~12 x 75~97	2,500	580
RA-1004	10ml x 4	15~16.5 x 104~110	2,500	650