

無刷數位離心機 Brushless Digital Type



CB-2200

特點：

- 採無碳刷馬達，無污染、防干擾。
- 數位式轉速及時間顯示。
- 10組記憶系統，可預先設定常用參數。
- ABS外殼及觀測視窗。
- 多重安全防護功能：不平衡斷電裝置、門蓋開啟具斷電裝置、自動安全門扣...等。

Feature :

- Brushless motor. No pollution and no interference
- Digital speed and time display
- 10 memory programs for presetting parameters
- ABS hard shell and see-through window
- Multiple safety devices included auto power off device when the rotor is imbalanced or door is opened during operation, auto-door-lock, etc

Model	CB-2200
轉速 Speed (rpm)	1000~8000
最大離心力 Max. RCF (xg)	6,660
最大容量 Max. Capacity (ml)	300
定時器 Timer (min:sec)	99:59
重量 Weight (kg)	8
尺寸 Size (W x L x H mm)	290 x 320 x 250
電源 Power	AC 110/220V, 50/60Hz

◎ 離心試管另購。 Tubes not included.

■ Rotor sets for CB-2200



Rotor	Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
RA-1508-1	15ml x 8	15~17.3 x 92~122	8,000	6,660
RA-1512S-1	15ml x 12	15~17.3 x 92~122	7,000	5,320
RA-5004	50ml x 4	27.5~29.5 x 84~117	7,500	6,040
RA-5006	50ml x 6	27.5~29.5 x 84~117	7,000	5,260

無刷數位離心機 Brushless Digital Type

20

離心機系列
Centrifuge Series



CB-2060

特點：

- 採無碳刷馬達，無污染、防干擾。
- 數位式轉速及時間顯示。
- 10組記憶系統，可預先設定常用參數。
- ABS外殼及觀測視窗。
- 多重安全防護功能：不平衡斷電裝置、門蓋開啟具斷電裝置、自動安全門扣...等。

Feature：

- Brushless motor. No pollution and no interference.
- Digital speed and time display
- 10 memory programs for presetting parameters
- ABS hard shell and see-through window
- Multiple safety devices included auto power off device when the rotor is imbalanced or door is opened during operation, auto-door-lock, etc

Model	CB-2060
轉速 Speed (rpm)	500~6000
最大離心力 Max. RCF (xg)	3,910
最大容量 Max. Capacity (ml)	300
定時器 Timer (min:sec)	99:59
重量 Weight (kg)	8
尺寸 Size (W x L x H mm)	290 x 320 x 250
電源 Power	AC 110/220V, 50/60Hz

◎ 離心試管另購。 Tubes not included.

■ Rotor sets for CB-2060



RA-1508



RA-1512S



RA-5004



RA-5006

Rotor	Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
RA-1508	15ml x 8	15~17.5 x 92~122	6,000	3,750
RA-1512S	15ml x 12	15~17.3 x 92~122	6,000	3,910
RA-5004	50ml x 4	27.5~29.5 x 84~117	6,000	3,870
RA-5006	50ml x 6	27.5~29.5 x 84~117	6,000	3,870