

中量型冷凍離心機 Refrigerated Type (Medium Capacity)

7

離心機系列 Centrifuge Series



特點：

- 採用無刷馬達。
- LCD數位面板，具溫度、轉速及離心力雙重切換、設定及顯示功能。
- 10組記憶系統，可預先設定常用參數。
- 不平衡斷電、過載、故障指示安全防護裝置。
- 具預冷功能。

Feature：

- Brushless Motor
- LCD digital panel with temperature, RPM & RCF switch, setting and display functions
- 10 memory programs for pre-setting parameters
- With auto power-off safety devices included rotor imbalance, overload, and error indication
- With the pre-cooling Function



掀蓋內部 Inside


Model	CN-350R	CN-350RD	CN-350RDS
轉速 Speed (rpm)	500~6,000		
最大離心力 Max. RCF (xg)	4,270		
最大容量 Max. Capacity (ml)	400		
溫度範圍 Temperature range	0°C~40°C		
定時器 Timer (min:sec)	99 : 59		
重量 Weight (kg)	85	55	60
尺寸 Size (W x L x H mm)	430 x 580 x 850	650 x 550 x 350	440 x 720 x 350
電源 Power	AC 110/220V, 50/60Hz		

◎ 離心試管另購。 Tubes not included.

■ Rotor sets for CN-350R, CN-350RD, CN-350RDS



Rotor	Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
RS-50	1.5ml x 48	9.5~10.5 x 38~50	3,500	1,320
	6ml x 36	11~12 x 85~106	3,500	1,940
	6ml x 40	11~12.5 x 63~87	3,500	1,630
	10ml x 24	14~15.5 x 87~108	3,500	1,980
	15ml x 8	15~17.3 x 90~115	3,500	1,980
	15ml x 12 (plastic)	15~17 x 89~115	3,500	2,000
	15ml x 16	15~17.3 x 90~110	3,500	1,980
	50ml x 4 (plastic)	26.5~29.5 x 90~120	4,000	2,560
	50ml x 4 (glass)	32~35 x 90~112	4,000	2,490
	100ml x 4	42~45 x 97~112	3,500	1,980



Rotor	Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
RS-55	1.5ml x 48	9.5~10.5 x 38~50	3,500	1,320
	6ml x 36	11~12 x 85~106	3,500	1,940
	6ml x 40	11~12.5 x 63~87	3,500	1,630
	10ml x 24	14~15.5 x 87~108	3,500	1,980
	15ml x 8	15~17.3 x 90~115	4,500	3,270
	15ml x 12 (plastic)	15~17 x 89~120	3,500	2,000
	15ml x 16	15~17.3 x 90~120	4,500	3,270
	50ml x 4 (plastic)	26.5~29.5 x 90~120	5,000	4,000
	50ml x 4 (glass)	32~35 x 90~112	5,000	3,890
	100ml x 4	42~45 x 97~112	4,500	3,270

■ Rotor sets for CN-350R, CN-350RD, CN-350RDS



Rotor	Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
RA-1024	10ml x 24	15~17.3 x 90~107	3,500	1,760
RA-1508-1	15ml x 8	15~17.3 x 92~130	6,000	3,750
RA-1512H-1	15ml x 12	15~17.3 x 92~122	6,000	4,270
RA-1516	15ml x 16	15~17.3 x 90~115	4,000	2,220



Rotor	Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
RA-5004	50ml x 4	27.5~29.5 x 84~135	6,000	3,870
RA-5006	50ml x 6	27.5~29.5 x 84~135	6,000	3,870
RA-5008	50ml x 8	27.5~29.5 x 84~116	6,000	4,270