

## 中量型離心機 Medium Capacity Type

15

離心機系列 Centrifuge Series



CN-350

### 特點：

- 採用無刷馬達。
- LCD數位面板，具轉速、離心力之雙重切換、設定及顯示功能。
- 10組記憶系統，可預先設定常用參數。
- 不平衡斷電、過載、故障指示安全防護裝置。

### Feature：

- Brushless Motor
- LCD digital panel with RPM & RCF switch, setting & display functions
- 10 memory programs for presetting parameters
- With auto power-off safety devices included rotor imbalance, overload, and error indication

| Model                   | CN-350               |
|-------------------------|----------------------|
| 轉速 Speed (rpm)          | 500~6,000            |
| 最大離心力 Max. RCF (xg)     | 4,270                |
| 最大容量 Max. Capacity (ml) | 400                  |
| 定時器 Timer (min:sec)     | 99:59                |
| 重量 Weight (kg)          | 20.5                 |
| 尺寸 Size (W x L x H mm)  | 360 x 480 x 340      |
| 電源 Power                | AC 110/220V, 50/60Hz |

◎ 離心試管另購。 Tubes not included.



CN-3302

### 特點：

- 採用無刷馬達。
- 數位式轉速顯示。
- 不平衡斷電、過載、故障指示安全防護裝置。
- 轉速及離心力設定功能。

### Feature：

- Brushless Motor
- Digital speed display
- With auto power-off safety devices included rotor imbalance, overload, and error indication
- Speed and RCF settings

| Model                   | CN-3302              |
|-------------------------|----------------------|
| 轉速 Speed (rpm)          | 500~6,000            |
| 最大離心力 Max. RCF (xg)     | 4,270                |
| 最大容量 Max. Capacity (ml) | 400                  |
| 定時器 Timer (min:sec)     | 99:59                |
| 重量 Weight (kg)          | 18                   |
| 尺寸 Size (W x L x H mm)  | 360 x 480 x 340      |
| 電源 Power                | AC 110/220V, 50/60Hz |

◎ 離心試管另購。 Tubes not included.

## ■ Rotor sets for CN-350, CN-3302



1.5ml x 48pcs



6ml x 36pcs




6ml x 40pcs



10ml x 24pcs



15ml x 8pcs

15ml x 12pcs  
(plastic)

15ml x 16pcs

50ml x 4pcs  
(plastic)50ml x 4pcs  
(glass)

100ml x 4pcs

| Rotor | Capacity (ml x pcs) | Tube Size (Ø x L) mm | Max. Speed (rpm) | Max. RCF (xg) |
|-------|---------------------|----------------------|------------------|---------------|
| RS-50 | 1.5ml x 48          | 9.5~10.5 x 38~50     | 3,500            | 1,320         |
|       | 6ml x 36            | 11~12 x 85~106       | 3,500            | 1,940         |
|       | 6ml x 40            | 11~12.5 x 63~87      | 3,500            | 1,630         |
|       | 10ml x 24           | 14~15.5 x 87~108     | 3,500            | 1,980         |
|       | 15ml x 8            | 15~17.3 x 90~115     | 3,500            | 1,980         |
|       | 15ml x 12 (plastic) | 15~17 x 89~115       | 3,500            | 2,000         |
|       | 15ml x 16           | 14~16.6 x 90~110     | 3,500            | 1,980         |
|       | 50ml x 4 (plastic)  | 26.5~29.5 x 90~120   | 4,000            | 2,560         |
|       | 50ml x 4 (glass)    | 32~35 x 90~112       | 4,000            | 2,490         |
|       | 100ml x 4           | 42~45 x 97~112       | 3,500            | 1,980         |



1.5ml x 48pcs



6ml x 36pcs



6ml x 40pcs



10ml x 24pcs



15ml x 8pcs

15ml x 12pcs  
(plastic)

15ml x 16pcs

50ml x 4pcs  
(plastic)50ml x 4pcs  
(glass)

100ml x 4pcs

| Rotor | Capacity (ml x pcs) | Tube Size (Ø x L) mm | Max. Speed (rpm) | Max. RCF (xg) |
|-------|---------------------|----------------------|------------------|---------------|
| RS-55 | 1.5ml x 48          | 9.5~10.5 x 38~50     | 3,500            | 1,320         |
|       | 6ml x 36            | 11~12 x 85~106       | 3,500            | 1,940         |
|       | 6ml x 40            | 11~12.5 x 63~87      | 3,500            | 1,630         |
|       | 10ml x 24           | 14~15.5 x 87~108     | 3,500            | 1,980         |
|       | 15ml x 8            | 15~17.3 x 90~115     | 4,500            | 3,270         |
|       | 15ml x 12 (plastic) | 15~17 x 89~120       | 3,500            | 2,000         |
|       | 15ml x 16           | 14~16.6 x 90~120     | 4,500            | 3,270         |
|       | 50ml x 4 (plastic)  | 26.5~29.5 x 90~120   | 5,000            | 4,000         |
|       | 50ml x 4 (glass)    | 32~35 x 90~112       | 5,000            | 3,890         |
|       | 100ml x 4           | 42~45 x 97~112       | 4,500            | 3,270         |

## ■ Rotor sets for CN-350, CN-3302



RA-1024



RA-1508



RA-1512H



RA-1516

| Rotor    | Capacity (ml x pcs) | Tube Size (Ø x L) mm | Max. Speed (rpm) | Max. RCF (xg) |
|----------|---------------------|----------------------|------------------|---------------|
| RA-1024  | 10ml x 24           | 15~17.3 x 90~107     | 3,500            | 1,760         |
| RA-1508  | 15ml x 8            | 14~16.6 x 92~130     | 6,000            | 3,750         |
| RA-1512H | 15ml x 12           | 14~16.6 x 92~122     | 6,000            | 4,270         |
| RA-1516  | 15ml x 16           | 15~17.3 x 90~112     | 4,000            | 2,220         |



RA-5004



RA-5006



RA-5008

| Rotor   | Capacity (ml x pcs) | Tube Size (Ø x L) mm | Max. Speed (rpm) | Max. RCF (xg) |
|---------|---------------------|----------------------|------------------|---------------|
| RA-5004 | 50ml x 4            | 27.5~29.5 x 84~135   | 6,000            | 3,870         |
| RA-5006 | 50ml x 6            | 27.5~29.5 x 84~135   | 6,000            | 3,870         |
| RA-5008 | 50ml x 8            | 27.5~29.5 x 84~116   | 6,000            | 4,270         |