

簡易型離心機 Simple Type

23

離心機系列 Centrifuge Series



CN-810



CN-820 / CN-830

特點：

- 採用無刷馬達。
- ABS外殼及觀測視窗。
- 可更換6ml用試管套。

Feature :

- Brushless Motor
- ABS hard shell and see-through window
- Also available to use the adapters of 6ml

Model	CN-810	CN-820	CN-830
轉速 Speed (rpm)	Fixed 3,000(50Hz) / 3,500(60Hz)	0~3,000 (50Hz) / 0~3,500 (60Hz)	
最大離心力 Max. RCF (xg)	1,040 / 1,420		
最大容量 Max. Capacity (ml)	10 x 8		10 x 12
定時器 Timer	0~30min or “ON”		
重量 Weight (kg)	8		
尺寸 Size (W x L x H mm)	290 x 320 x 250		
電源 Power	AC 110V / 220V, 50/60Hz		

◎ 離心試管另購。 Tubes not included.

■ Rotor sets for CN-810, CN-820, CN-830



R-8



R-12

Rotor		Capacity (ml x pcs)	Tube Size (Ø x L) mm	Max. Speed (rpm)	Max. RCF (xg)
CN-810 CN-820	R-8	6ml x 8	10~12.5 x 75~107	3,500	1,270
		10ml x 8	12~16.5 x 95~117		1,420
CN-830	R-12	6ml x 12	10~12.5 x 75~90		1,270
		10ml x 12	12~16.5 x 88~100		1,420