

# 溶出試驗器 Dissolution Tester

100

其他泛用儀器、設備

Other Universal Instrument & Equipment



DT-6

## 特點：

- 本機可用於試驗藥品於定溫下之溶解程度，其轉速採無段變速方式，可調整速度，數位式顯示轉速。反應水槽內的溫度可依需要設定，並採自動循環裝置，以維持水槽溫度均等。驅動部採電動式上下昇降，自動定位停止。加裝防護裝置，及液位開關裝置，以確保使用安全。攪拌棒採SUS316製造。

## Feature：

- It is used to test the dissolution condition of chemicals under the set temperature. It's rotary speed is non-limetary variation and can be adjusted at any required speed. The digital readout provides a continuously updated accurate reading of shaft rpm. The temperature of water-baths can be set to a need one, in special, it can keep up the set temperature automatically. The circuit breaker device can prevent from leakage current. Thus, it can provide safety when use. The paddle shafts are made of SUS 316.

Model	DT-6
懸籃 Basket	6
轉速 Speed (rpm)	數位轉速計 Digital Tachometer 10 ~ 220
加熱管 Pipe Heater	0.9kw
尺寸 Dimension (mm)	935 x 350 x 640
重量 Weight (kg)	62
電源 Power	AC 220V, 50/60Hz

## 配件 / Accessories：



圓底燒杯  
Round Bottom Flask  
6 PCS



杯蓋  
Cover  
6 PCS



攪拌棒 / 槳式攪拌棒  
Rotation Shaft  
Paddle Shaft  
6 / 6 PCS



線籃  
Sinker  
6 PCS



網籃  
Basket  
6 PCS