


## 磨擦式粉碎機 Frictional Efficient Mill



FTM-30



刀片內部  
A close-up photo  
of the frictional  
disc set.

刀具外部  
A close-up photo  
of the thickness  
control bars.

### 特點：

- 本機採用特殊鑄鋼製成。
- 本機利用互相摩擦，使試料達到粉碎。
- 粉碎細度可任意調節。
- 連續式粉碎。
- 適用於礦石、煤炭、焦炭等粉碎。

### Feature :

- The machine body is made of cast steel.
- By the mutual friction to mill materials.
- Particle fineness is adjustable.
- Continuous output during operation.
- Applies to mill mineral · coal carbon... · etc.

Model	FTM-30
粉碎細度 Size of particles	≤ 0.3m/m
粉碎量 Capacity	3 kg / H
入料口尺寸 Inlet Size (mm)	W100 x D70 mm
出料口尺寸 Outlet Size (mm)	W45 x D25 mm
外部尺寸 Size (mm)	W360 x D410 x H870
電源 Power	AC110 / 220V, 50 / 60Hz

## 自動乳鉢破碎機 Auto-Pestle Crusher



AWD-200

### 特點：

- 乳棒以自由落體的方式錘擊物料，進行破碎、研磨。
- 錘擊的速度可調整。
- 數位式錘擊累計顯示，並可設定錘擊次數。
- 乳鉢與乳棒為不鏽鋼材質。

### Feature :

- Pestle in free fall way to beat materials in the mortar.
- Beating speed is adjustable.
- Digital display includes the cumulative.
- The mortar and pestle are stainless steel material.

Model	AWD-200
馬力 Motor	90W
速度 Speed (epm)(50/60Hz)	0~120(50Hz)/0~100(60Hz)
行程 Stroke(mm)	60
設定次數 Beating Frequency	1~999,999(6位數顯示)
重量 Weight (kg)	39
尺寸 Size (W x L x H mm)	385 x 310 x 475
電源 Power	AC100V/110V/220V, 50/60Hz