

電磁攪拌機 Magnetic Stirrer

70

電磁攪拌機系列
Magnetic Stirrer Series



MA-30



MA-626



MS-203

| Model | MA-30 | MA-626 | MS-203 |
|---|-----------------------|-------------------|----------------|
| Stirring Capacity (ml) (H ₂ O) | 100~5,000 | 100~2,000 x 6 | 50~3,000 |
| Speed (rpm) | 700~1,500 | | 90~1,500 (LED) |
| Plate Size (mm) | 180 x 180 | 300 x 450 | Ø150 |
| Plate Material | SUS 304 Stainless | SUS 304 Stainless | Aluminum |
| Power | AC 110V/220V, 50/60Hz | | |



JTR-500



JS-M



JS-B

| Model | JTR-100 | JTR-300 | JTR-500 |
|------------------------|-----------------|-----------------|-----------------|
| Stirring Capacity (ml) | 1,000 | 3,000 | 3,000 |
| Speed (rpm) | 100~1,500 | 100~1,500 | 100~1,500 |
| Plate Size (mm) | 110 x 110 | 160 x 160 | 220 x 220 |
| Plate Material | ABS | ABS | ABS |
| Accessory | Stirrer bar x 1 | Stirrer bar x 1 | Stirrer bar x 1 |

| Model | JS-M | JS-B |
|-------------------------|----------------|-----------------|
| Stirring Capacity (ml) | 50~2,000 | 50~3,000 |
| Speed (rpm) | Max. 1,500 | Max. 1,500 |
| Plate Size (mm) | 110 x 100 | 200 x 200 |
| Body Size (W x L x Hmm) | 165 x 105 x 90 | 285 x 200 x 125 |

Power : AC 110V , 60Hz