



Specifications:

MODEL	P3-50	P20-300	TP-300	P-230/3500
Display	digital setting display		LCD display with touch screen	digital setting display
Working range (rpm)	3 ~ 50	20 ~ 300	20 ~ 300	20 ~ 230
Flow rate (ml/min)	As shown below			
Dimension (cm)	W21 / D31 / H14	W17 / D26 / H23	W26 / D33 / H17	
Power source	110/220V	60/50Hz	1A	110/220V 60/50Hz 2A

Special function

1. 具自動流量校正功能 (TP-300)
2. 具定量定時饋料功能
3. 可定時定量正逆轉功能
4. 可擴充 4-20mA 由外部類比信號控制
5. 可擴充 RS-485 由外部數位信號控制

1.6mm Wall Thickness Tubing							
bore \ Flow rate (ml/min)	0.8mm	1.6mm	3.2mm	4.8mm	6.4mm	8.0mm	9.6mm
P3-50	0.21~3.5	0.81~13.5	3 ~ 50	6.6 ~110	10.8 ~ 180	15 ~ 250	—
P20-300 & TP-300	1.4 ~ 21	5.4 ~ 81	20 ~ 300	44 ~ 660	72 ~ 1080	100 ~ 1500	—
P-230/3500	2 ~ 24	8 ~ 87	36 ~ 390	80 ~ 870	140 ~ 1500	220 ~ 2400	320 ~ 3500