

# 大型振盪培養機 Giant Shaking Incubator



TRL-1000 / 1200

## 特點：

- 可用於一般實驗、振盪、生物科技、遺傳工程、動植物細胞及微生物培養用。本機各支軸偏心迴轉，有效的拋物線方式振盪瓶內試料，五點式平衡振盪運轉安靜順暢，可設置恆溫恆濕培養室。

## Feature :

- Suitable to laboratory, shaking, biological tech, genetic, engineering, animal and plant cell culture. Each support axis rotates centripetally which can effectively shake vials materials by parabolic. 5-point shaking evenly and quietly, can be installed in constant temperature and humidity room.

Model	TRL-1000	TRL-1200
層板尺寸 Plate Size	(L)1000 x (W)1000mm	(L)1200 x (W)1200mm
尺寸 Dimension	(L)1085 x (W)1205mm	(L)1285 x (W)1405mm
轉速 Speed	0 ~ 250rpm	
行程 Stroke	70mm	
重量 Weight	約 800 Kg	
電源 Power	3Φ AC220V, 50/60Hz	

Model	TRL-1000 Clamp Quantity								
Flask	50ml	125ml	250ml	300ml	500ml	1000ml	2000ml	3000ml	5000ml
1-Layer-Type	256	144	110	100	64	36	28	20	16
2-Layer-Type	256	144	110	100	64	36	28	20	16
3-Layer-Type	256	144	110	100	64	--	--	--	--

Model	TRL-1200 Clamp Quantity								
Flask	50ml	125ml	250ml	300ml	500ml	1000ml	2000ml	3000ml	5000ml
1-Layer-Type	342	210	156	144	100	64	36	28	20
2-Layer-Type	342	210	156	144	100	64	36	28	20
3-Layer-Type	342	210	156	144	100	--	--	--	--

# 分液漏斗振盪機 Funnel Shaker

88

振盪機系列  
Shaker Series



VD-12



VD-12-2S



V-U



VS-6

## 特點：

- 本機是將多種試料在同一條件下振盪攪拌均勻，適用於有機與無機化學、醫學、藥品等之基礎研究，或者臨床之溶媒抽出培養產業、排水的水質檢查等各種試驗。
- 指針式振盪次數指示錶。
- 旋鈕式定時器0~60 Min或數位式定時器0~99Hr, 59 Min 可供選擇。

## Feature :

- This shaker is efficient in shaking reagents of a variety of materials under same control parameters. It's suitable for organic/inorganic chemical, medical, pharmaceutical research, extraction and development of media, and waste water analysis.
- An analog meter for EPM indication.
- A mechanical timer (0~60Min) or digital timer (0~99Hr,59Min) is for option.

Model	VS-6	VD-12	V-U	VD-12S-2S
Oscillations (epm)	0 ~ 350		0 ~ 350 (Vertical) 0 ~ 300 (Horizontal)	0 ~ 300
Capacity & Funnel Holder Q'ty	100 x 6	100 x 12	100 x 6	100 x 12
	200 x 6	200 x 12	200 x 6	200 x 12
	300 x 4	300 x 8	300 x 4	300 x 8
	500 x 4	500 x 8	500 x 4	500 x 8
	1000 x 3	1000 x 6	1000 x 3	1000 x 6
	---	---	---	2000 x 4
Size (W x L x H mm)	430 x 430 x 560	480 x 430 x 560	455 x 430 x 560	430 x 710 x 730
Weight (Kg)	59 Kg	65 Kg	71 Kg	81 Kg

Power : AC 110V/220V, 50/60Hz

## 分液漏斗振盪機配件 Funnel Shaker Accessories

配件(選購品) / Accessories : (For optional purchases)



適用之容器 / Applicable vessels

