

彈壓式升降高速攪拌機 Spring-lift & high-speed stirrer

46

攪拌機系列
Stirrer Series



ST-202



ST-202

特點：

- 數位式轉速顯示。
- 彈壓式升降，操作迅速輕便。
- 攪拌棒可上下調整。
- 無段變速，轉速穩定，具扭力補償功能。
- 過電流保護、過熱保護及保險絲安全裝置。
- 精密結構，運轉安靜。
- 攪拌棒、夾頭及軸心皆採不鏽鋼製。

Feature：

- Digital speed (rpm) display
- Spring-lift pole for easy operation
- Stirring shaft is adjustable up and down
- Stepless and stable speed; auto torque-compensated function
- With the over-current, over-heat and fuse protection devices
- Precise structure and quiet operation
- Stainless stirring shaft、chuck and axis

Model	ST-202
馬力 Motor	DC 210W
轉速 Speed (rpm)	200~3,000
扭力 Torque (kgf·cm / N·m)	7 / 0.69
重量 Weight (kg)	10
機身尺寸 Body Size (W x L x H mm)	130 x 283 x 273
整機尺寸 Overall dimension (W x L x H mm)	330 x 303 x 695
電源 Power	AC 110V/220V, 50/60Hz

- 標準配件
- 攪拌棒：Ø8 x 500mm
- 攪拌葉：No.11 Ø50mm
- 腳座：W328 x L262 x H55mm · 3.3kg
- 選購配件 (詳見P50~P52)

- Standard accessories
- Stirring shaft：Ø8 x 500mm
- Propeller：No.11 Ø50mm
- Stand base：W328 x L262 x H55mm · 3.3kg
- Optional purchasing accessories. (Ref. P50~P52)